

# Topic 7. Molecular weight measurement

Expected learning outcomes:

By the end of this lesson, student should be able to...

- discuss and explain the experimental method to measure the molecular weight distribution of polymers
- discuss and explain basic principles of GPC operation
- understand and discuss the rationale for separation of polymer chains by size and the role of detectors in GPC
- apply methods of calibration to GPC, including the universal calibration curve

## **Topic 7. Molecular weight measurement**

### **Sub-Topics:**

7.1 Principles of gel permeation chromatography (GPC)

7.2 Parts of GPC and calibration



7.3 Universal calibration curve

7.3.1 Universal calibration curve – calibration with standards of same polymer

7.3.2 Universal calibration curve – calibration with standards of different polymer


7.3.3 Universal calibration curve – calibration with viscosity

## Separation based on molecular weight




- The smaller polymer chains take longer time than the larger polymer chains to go through all the pores.
- The polymer chains with small size elute at longer times.

## Gel Permeation Chromatography (GPC):




- The smaller polymer chains take longer time than the larger polymer chains to go through all the pores.
- The polymer chains with small size elute at longer times (or larger elution volumes).


## 7.2 Parts of GPC and Calibration




## Gel Permeation Chromatography (GPC) - Calibration:



## Gel Permeation Chromatography (GPC) - Calibration:



## Gel Permeation Chromatography (GPC):





Detector: Differential Refractive Index (DRI)

DRI detector depends on the mass concentration:

$$C_m = \frac{N \times M_P}{V_{cell}} \quad (V_{cell} = \text{volume of the DRI detection cell})$$


$V_{cell} = 8 \mu\text{L}$ , very small, so  $C_m = N \times M_p$

### 7.3.1 Universal calibration curve – calibration with standards of same polymer


#### Gel Permeation Chromatography (GPC) – $M_n$ :

$$Int\#1 = \sum \frac{w_i}{M_i} = \sum \frac{\frac{N_i M_i}{M_i}}{\sum \frac{N_i M_i}{M_i}} = \frac{1}{\sum N_i M_i} \sum \frac{N_i M_i}{M_i} = \frac{\sum N_i}{\sum N_i M_i}$$

$$Int\#1 = M_n^{-1}$$



$$\ln(M_p) = m \times V_e + p$$





$$\frac{I_{DRI}}{\sum I_{DRI}} = \frac{N_i M_i / V_{cell}}{\sum N_i M_i / V_{cell}} = \frac{N_i M_i}{\sum N_i M_i} = w_i$$


## Gel Permeation Chromatography (GPC) – $M_w$ and PDI:

$$Int\#1 = \sum \frac{w_i}{M_i} = \sum \frac{\frac{N_i M_i}{\sum N_i M_i}}{M_i} = \frac{1}{\sum N_i M_i} \sum \frac{N_i M_i}{M_i} = \frac{\sum N_i}{\sum N_i M_i}$$

$$Int\#1 = M_n^{-1}$$




$$\ln(M_p) = m \times V_e + p$$



$$Int\#2 = \sum w_i M_i = \sum \frac{N_i M_i}{\sum N_i M_i} \times M_i = \frac{\sum N_i M_i^2}{\sum N_i M_i}$$

$$Int\#2 = M_w$$

$$Int\#1 \times Int\#2 = M_w / M_n = PDI$$




## Gel Permeation Chromatography (GPC) – Limitations:

- ⇒ The calibration curve is only valid for the same type of polymer used for the calibration, under the same conditions (solvent, temperature, column).
- ⇒ Narrow molecular weight distributions are only available for a limited number of polymers. Nevertheless, “apparent” molecular weights based on a polystyrene calibration curve are routinely reported in the literature.

### 7.3.2 Universal calibration curve – calibration with standards of different polymer

Two different polymers with the same hydrodynamic volume elute at the same elution volume.



## Gel Permeation Chromatography (GPC) – Universal Calibration:

$$[\eta] = 2.5 \times N_A \times \frac{V_h}{M_P}$$

$$[\eta] \times M_P = 2.5 \times N_A \times V_h$$

According to the Mark-Houwink-Sakurada equation:

$$[\eta] = K \times M_P^a \text{ with } 0.5 < a < 0.8$$

$$[\eta] \times M_P = K \times M_P^{a+1} = 2.5 \times N_A \times V_h$$

## Gel Permeation Chromatography (GPC) – Universal Calibration:

$$[\eta] \times M_P = K \times M_P^{a+1} = 2.5 \times N_A \times V_h$$

Two different polymers *A* and *B* eluting with the same elution volume have the same hydrodynamic volume  $V_h$  so that:

$$K_A M_{PA}^{a_A+1} = 2.5 \times N_A \times V_h (A \text{ or } B) = K_B M_{PB}^{a_B+1}$$

**A, B with  
known  $\alpha$  and K**


$$M_{PB} = \left( \frac{K_A}{K_B} \right)^{\frac{1}{a_B+1}} M_{PA}^{\frac{a_A+1}{a_B+1}}$$

**constants**

$$\ln(M_{PB}) = \frac{1}{a_B+1} \ln\left(\frac{K_A}{K_B}\right) + \frac{a_A+1}{a_B+1} \ln(M_{PA})$$


## Gel Permeation Chromatography (GPC) – Universal Calibration:

**constants**


$$\ln(M_{PB}) = \frac{1}{a_B + 1} \ln\left(\frac{K_A}{K_B}\right) + \frac{a_A + 1}{a_B + 1} \ln(M_{PA})$$

If the column has been calibrated with polymer A, then we have:

$$\ln(M_{PA}) = m \times V_e + p$$


$$\ln(M_{PB}) = \left[ \frac{1}{a_B + 1} \ln\left(\frac{K_A}{K_B}\right) + \frac{a_B + 1}{a_A + 1} \times p \right] + \frac{a_A + 1}{a_B + 1} \times m \times V_e$$

**constant =  $p'$**

**constant =  $m'$**



$$\ln(M_{PB}) = m' \times V_e + p'$$




## Gel Permeation Chromatography (GPC) – $M_n$ , $M_w$ , and PDI:

$$Int\#1 = \sum \frac{w_i}{M_i} = \sum \frac{\frac{N_i M_i}{\sum N_i M_i}}{M_i} = \frac{1}{\sum N_i M_i} \sum \frac{N_i M_i}{M_i} = \frac{\sum N_i}{\sum N_i M_i}$$

$$Int\#1 = M_n^{-1}$$




$$\ln(M_{PA}) = m \times V_e + p$$



$$Int\#2 = \sum w_i M_i = \sum \frac{N_i M_i}{\sum N_i M_i} \times M_i = \frac{\sum N_i M_i^2}{\sum N_i M_i}$$

$$Int\#2 = M_w$$



$$Int\#1 \times Int\#2 = M_w / M_n = PDI$$




# Gel Permeation Chromatography (GPC) – $M_n$ , $M_w$ , and PDI:

$$Int\#1 = \sum \frac{w_i}{M_i} = \sum \frac{\frac{N_i M_i}{\sum N_i M_i}}{M_i} = \frac{1}{\sum N_i M_i} \sum \frac{N_i M_i}{M_i} = \frac{\sum N_i}{\sum N_i M_i}$$

$$Int\#1 = M_n^{-1}$$




$$\ln(M_{PB}) = m' \times V_e + p'$$



$$Int\#2 = \sum w_i M_i = \sum \frac{N_i M_i}{\sum N_i M_i} \times M_i = \frac{\sum N_i M_i^2}{\sum N_i M_i}$$

$$Int\#2 = M_w$$

$$Int\#1 \times Int\#2 = M_w/M_n = PDI$$





### 7.3.3 Universal calibration curve – calibration with viscosity

If the  $\alpha$  and  $K$  for B is unknown,  $m'$  and  $p'$  are not available for the calibration

$$\ln(M_{PB}) = m' \times V_e + p'$$

Viscosity detector



The concentration of the sample in the cell is very low

$$\text{So } [\eta] = \frac{1}{C_m} \frac{\eta - \eta_o}{\eta_o}$$

## Gel Permeation Chromatography (GPC) – Universal Calibration:

$$[\eta] \times M_P = K \times M_P^{a+1} = 2.5 \times N_A \times V_h$$


Two different polymers  $A$  and  $B$  eluting with the same elution volume have the same hydrodynamic volume  $V_h$  so that:

$$[\eta]_B M_{PB} = 2.5 \times N_A \times V_h (A \text{ or } B) = K_A M_{PA}^{a_A+1}$$

experimentally known quantities

$$M_{PB} = \frac{K_A M_{PA}^{a_A+1}}{[\eta]_B}$$


## Gel Permeation Chromatography (GPC) – Universal Calibration:





## Gel Permeation Chromatography (GPC) – $M_n$ , $M_w$ , and PDI:

$$Int\#1 = \sum \frac{w_i}{M_i} = \sum \frac{\frac{N_i M_i}{\sum N_i M_i}}{M_i} = \frac{1}{\sum N_i M_i} \sum \frac{N_i M_i}{M_i} = \frac{\sum N_i}{\sum N_i M_i}$$

$$Int\#1 = M_n^{-1}$$




$$M_{Bi} = \frac{K_A M_{PA}^{a_A+1}}{[\eta]_{Bi}}$$



$$Int\#2 = \sum w_i M_i = \sum \frac{N_i M_i}{\sum N_i M_i} \times M_i = \frac{\sum N_i M_i^2}{\sum N_i M_i}$$


$$Int\#2 = M_w$$

$$Int\#1 \times Int\#2 = M_w / M_n = PDI$$




# Gel Permeation Chromatography (GPC) – Universal Calibration:


Benoit et al. *Polym. Lett.* **1967**, 5, 753




linear



comb



star



The molecular weight of a new polymer can be measured by GPC using universal calibration; but the sample has to be prepared with an accurate concentration.

# Summary / Review

- Gel permeation chromatography separates chains by hydrodynamic volume
  - Only method to produce entire molecular weight distribution (MWD) curve
- Calibration curve correlates elution volume and molecular weight peak ( $M_p$ )
- Refractive index and viscometer most popular detectors
- Universal calibration curve allows measurement of MWD curve using calibration made with a different polymer (when  $\alpha$  and  $K$  are known)



## Summary of Methods:

Ebullioscopy and cryoscopy:  $M_n < 30,000 \text{ g.mol}^{-1}$ , yields  $M_n$ .

Vapor pressure depression:  $M_n < 20,000 \text{ g.mol}^{-1}$ , yields  $M_n$ .

Membrane osmometry: any  $M_n$ .

End group analysis:  $M_n < 1,000,000 \text{ g.mol}^{-1}$ , yields  $M_n$ ; depends on the efficiency of the labeling reaction.

Light scattering: any  $M_n$ , yields  $M_w$ .

Viscometry: any  $M_n$ , yields  $M_v \approx M_w$

Gel permeation chromatography:  $M_n < 10,000,000 \text{ g.mol}^{-1}$ , yields  $M_n$ ,  $M_w$ , PDI, and the entire distribution.